**kasios®**

# kasios® TCP

**Sustituto Oseo Sintético**  
**Synthetic Bone Substitute**

*Dental*

## Composición / Composition

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*Synthetic phosphocalcic ceramic made from  $\beta$  tri-calcium phosphate ( $\beta$ TCP) which closely resembles the mineral phase of natural bone.*

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*Kasios®TCP is entirely synthetic, so that risks as regards immune response or infection are eliminated.*

*Cavity filling and graft integration radiological monitoring is not impeded since the TCP granules are radio-opaque.*

*The shape of Kasios®TCP granules solves the problem of filling irregularly shaped cavities. Because of its composition, Kasios®TCP is biocompatible. Osteo-integration and resorption proceed at a rapid pace.*

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*Hold the blister firmly in one hand and slowly peel off the seal, being careful to hold at a slant to avoid spilling granules.*

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*Mix the Kasios<sup>®</sup> TCP granules with either the patient's blood, saline solution, or distilled water. Do not humidify the granules excessively. Any surplus liquid should be wiped away with the corner of a sterile sponge.*

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*The flat side of a sterile spatula may be used to mix Kasios<sup>®</sup> TCP granules and the liquid. Humidified granules can be transported to the surgical area using the curved part of the spatula.*

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*Gently fill each defect progressively up to the highest level of the bone cavity. Osteo-integration will be improved if blood and vessels can penetrate the graft.*

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*Once filling with the Kasios<sup>®</sup> TCP granules is completed, the flap should be carefully adjusted and completely closed. A membrane can also be used to ensure full closure.*

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